

IκB-alpha(N-term) Monoclonal Antibody

Catalog No: #27065

Package Size: #27065-1 50ul #27065-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	IκB-alpha(N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu P Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids of human NFKBIA
Conjugates	Unconjugated
Target Name	IκB-alpha
Other Names	NF-kappa-B inhibitor alpha <input type="checkbox"/> IκB-alpha <input type="checkbox"/> IκappaBalpha <input type="checkbox"/> kappa-B-alpha <input type="checkbox"/> nuclear factor of kappa light chain gene enhancer in B-cells <input type="checkbox"/> major histocompatibility complex enhancer-binding protein MAD3
Accession No.	UniProt:P25963 NCBI Gene ID:4792
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

Application Details

Predicted MW: 36kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected NFKBIA on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.