

TAB1(N-term) Monoclonal Antibody

Catalog No: #27075



Package Size: #27075-1 50ul #27075-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TAB1(N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu D Ms P Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids of human TAB1
Target Name	TAB1
Other Names	TGF-beta-activated kinase 1 and MAP3K7-binding protein 1 □ TAK1-binding protein 1 □ transforming growth factor beta-activated kinase-binding protein 1 □ mitogen-activated protein kinase kinase kinase 7-interacting protein 1
Accession No.	UniProt:Q15750 NCBI Gene ID:10454
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 55kd
Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected TAB1 on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.