

HNRNPU Monoclonal Antibody

Catalog No: #27076



Package Size: #27076-1 50ul #27076-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HNRNPU Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Rb B P
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids of human HNRNPU
Target Name	HNRNPU
Other Names	Scaffold attachment factor A
Accession No.	UniProt:Q00839 NCBI Gene ID:3192
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 90kd
Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected HNRNPU on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.