

COX4(isoform1) Monoclonal Antibody

Catalog No: #27080



Package Size: #27080-1 50ul #27080-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	COX4(isoform1) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC FC
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	A synthetic peptide (NPIQLGLASKWDYEKNEWKKC) of human COX4(C-terminus).
Target Name	COX4(isoform1)
Other Names	cytochrome c oxidase subunit 4 isoform 2, mitochondrial COX IV-2 cytochrome c oxidase subunit IV-like 2
Accession No.	UniProt:Q96KJ9 NCBI Gene ID:84701
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C


Application Details

Predicted MW: 20kd
Western blotting: 1:1000


Images



1:4000 dilution of this antibody detected COX (C-terminus)on 20



1:150 dilution from a previous lot detected COX4 (C-terminus) in Hela



K562 cells stained with COX4 (red, 1:100 dilution), followed by FITC-conjugated goat anti-mouse IgG. Blue line histogram represents the isotype control, normal mouse IgG

Note: This product is for in vitro research use only and is not intended for use in humans or animals.