

# LRRFIP1 Monoclonal Antibody

Catalog No: #27113



Package Size: #27113-1 50ul #27113-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                                             |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name          | LRRFIP1 Monoclonal Antibody                                                                                                 |
| Host Species          | Mouse                                                                                                                       |
| Clonality             | Monoclonal                                                                                                                  |
| Clone No.             | 4F3-D8-C9                                                                                                                   |
| Isotype               | IgG1                                                                                                                        |
| Purification          | Affinity purified                                                                                                           |
| Applications          | ELISA WB                                                                                                                    |
| Species Reactivity    | Hu                                                                                                                          |
| Immunogen Type        | Recombinant protein                                                                                                         |
| Immunogen Description | Purified recombinant human LRRFIP1 protein fragments expressed in E.coli.                                                   |
| Target Name           | LRRFIP1                                                                                                                     |
| Other Names           | GCF2, TRIP, FLAP1, GCF-2, FLAP-1, HUFI-1, FLIAP1                                                                            |
| Formulation           | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol |

## Application Details

Predicted MW: 160kd  
Western blotting: 1:1000

## Images



Western blot detection of LRRFIP1 in Raji,C2C12,MCF7&Jurkat cell lysates using LRRFIP1 antibody (1:1000 diluted).  
Predicted band size:89KDa.  
Observed band size:160KDa

## Background

Transcriptional repressor which preferentially binds to the GC-rich consensus sequence (5'-AGCCCCCGGCG-3') and may regulate expression of TNF, EGFR and PDGFA. May control smooth muscle cells proliferation following artery injury through PDGFA repression.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.