

FOXP3 Monoclonal Antibody

Catalog No: #27179



Package Size: #27179-1 50ul #27179-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FOXP3 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Affinity Purified
Applications	WB;IF/ICC
Species Reactivity	Human;Mouse;Rat
Specificity	Transfected Only
Immunogen Type	Recombinant Protein
Immunogen Description	Purified recombinant human FOXP3 protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	FOXP3
Other Names	AIID; AIID; DIETER; DIETER; Forkhead box P3; Forkhead box protein P3; foxp3; foxp3; FOXP3_HUMAN; FOXP3delta7; Immune dysregulation polyendocrinopathy enteropathy X linked; Immunodeficiency polyendocrinopathy enteropathy X nked; IPEX; IPEX; JM2; MGC141961
Accession No.	Uniprot: Q9BZS1 Gene ID: 50943
SDS-PAGE MW	47 kDa
Formulation	Ascites containing with 0.02% sodium azide and 50% glycerol.
Storage	store at -20 Λ C


Application Details

WB 1:500-1:1000; IF 1:50-1:200

Images



Immunocytochemistry staining of CHO-K1 cells which transfected by FOXP3 fixed in 1% Paraformaldehyde and using anti-FOXP3 antibody (dilution 1:100).



Western blot detection of FOXP3 antibody in CHO-K1 cell lysate. O A o and CHO-K1 transfected by FOXP3 o O B o cell lysate using FOXP3 antibody (1:300 diluted). Predicted band size: 47KDa. Observed band size: 50KDa.

Background

Probable transcription factor. Plays a critical role in the control of immune response.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.