

Phostensin Monoclonal Antibody

Catalog No: #27199



Package Size: #27199-1 50ul #27199-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Phostensin Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4B12-E11-D9
Isotype	IgG1
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of PPP1R18 and does not cross-react with related proteins.
Immunogen Type	Recombinant Protein
Immunogen Description	Purified recombinant human PPP1R18 protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	Phostensin
Other Names	HKMT1098; KIAA1949; Protein phosphatase 1 F actin cytoskeleton targeting subunit; PPP1R18;
Accession No.	Uniprot: Q6NYC8 Gene ID: 170954
SDS-PAGE MW	120kd
Formulation	Ascites
Storage	store at -20Λ C

Application Details

Western blotting: 1:200

Images



Western blot detection of PPP1R18 in HeLa,Raji,Ramos&Leucocyte cell lysates and using PPP1R18 antibody(1:200 diluted).Predicted band size: 120KDa
Observed band size: 120KDa.

Background

May target protein phosphatase 1 to F-actin cytoskeleton.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.