

ZSCAN4 Polyclonal Antibody

Catalog No: #27589



Package Size: #27589-1 50ul #27589-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZSCAN4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse
Immunogen Description	Recombinant fusion protein of human ZSCAN4 (NP_689890.1).
Conjugates	Unconjugated
Other Names	ZSCAN4; ZNF494; zinc finger and SCAN domain containing 4
Accession No.	Swiss-Prot#:Q8NAM6NCBI Gene ID:201516
Calculated MW	49kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZSCAN4 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.