

MICU2 Polyclonal Antibody

Catalog No: #27642



Package Size: #27642-1 50ul #27642-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MICU2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human MICU2 (NP_689939.1).
Other Names	MICU2; 1110008L20Rik; EFHA1; mitochondrial calcium uptake 2
Accession No.	Swiss-Prot#:Q8IYU8NCBI Gene ID:221154
Calculated MW	50kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MICU2 antibody.

Published Papers

el at., Quantitative proteomic analysis of scleras in guinea pig exposed to wavelength defocus. In J Proteomics on 2021 Jul 15 by Yingying Wen, Le Jin, et al..PMID:33964483, , (2021)

[PMID:33964483](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.