

AGTPBP1 Polyclonal Antibody

Catalog No: #27888

Package Size: #27888-1 50ul #27888-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	AGTPBP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human AGTPBP1 (NP_056054.2).
Conjugates	Unconjugated
Other Names	AGTPBP1; CCP1; NNA1; ATP/GTP binding protein 1
Accession No.	Swiss-Prot#:Q9UPW5NCBI Gene ID:23287
Calculated MW	138kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using AGTPBP1 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.