

SCD5 Polyclonal Antibody

Catalog No: #27911



Package Size: #27911-1 50ul #27911-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SCD5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human SCD5 (NP_001032671.2).
Other Names	SCD5; ACOD4; FADS4; HSCD5; SCD2; SCD4; stearyl-CoA desaturase 5
Accession No.	Swiss-Prot#:Q86SK9NCBI Gene ID:79966
Calculated MW	38kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using SCD5 antibody.

Background

Published Papers

el at., LncCCAT1 interaction protein PKM2 upregulates SREBP2 phosphorylation to promote osteosarcoma tumorigenesis by enhancing the Warburg effect and lipogenesis. In Int J Oncol on 2022 Apr by Feifei Pu, Jianxiang Liu,et al..PMID:35244192, , (2022)
[PMID:35244192](https://pubmed.ncbi.nlm.nih.gov/35244192/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.