

PID1 Polyclonal Antibody

Catalog No: #27951



Package Size: #27951-1 50ul #27951-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PID1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human PID1 (NP_060403.3).
Conjugates	Unconjugated
Other Names	PID1; HMFN2073; NYGGF4; P-CLI1; PCL11; phosphotyrosine interaction domain containing 1
Accession No.	Swiss-Prot#:Q7Z2X4NCBI Gene ID:55022
Calculated MW	27kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PID1 antibody.

Background

Published Papers

el at., Dual role of PID1 in regulating apoptosis induced by distinct anticancer-agents through AKT/Raf-1-dependent pathway in hepatocellular carcinomaInCell Death DiscovOn2023 Apr 28byJian Yang,Senlin Li et al..PMID:37117198, , (2023)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.