

# Rat Anti-mouse CD45/APC antibody

Catalog No: #28009



Package Size: #28009-1 50T #28009-2 100T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Anti-mouse CD45/APC antibody
Host Species	Rat
Clonality	Monoclonal
Purification	Affinity purified by Protein G
Applications	FC
Species Reactivity	Mouse
Conjugates	Unconjugated
Other Names	LCA; T200
Accession No.	Swiss-Prot#: P06800 NCBI Gene ID:19264
Formulation	0.01M PBS, 0.5%BSA, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.

## Application Details

Flow-Cyt=10ul/Test

Note: This product is for in vitro research use only and is not intended for use in humans or animals.