

Rat Anti-mouse CD206/FITC antibody

Catalog No: #28020

Package Size: #28020-1 50T #28020-2 100T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Anti-mouse CD206/FITC antibody
Host Species	Rat
Clonality	Monoclonal
Isotype	Rat IgG2a, k
Purification	Affinity purified by Protein G
Applications	FC
Species Reactivity	Mouse
Other Names	Mannose Receptor; C-type lectin domain family 13 member D; CD 206; CD206; CD206 antigen; CLEC13D; CLEC13DL; Macrophage mannose receptor; Mannose receptor C type 1; MMR; MRC 1; MRC-1; MRC1_MOUSE; Macrophage mannose receptor 1; MRC1L1; C-type lectin domain family 13 member D-like; Macrophage mannose receptor 1-like protein 1.
Accession No.	Swiss-Prot#: Q61830 NCBI Gene ID:17533
Formulation	0.01M PBS, 0.5%BSA, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.

Application Details

Flow-Cyt=10ul/Test

Note: This product is for in vitro research use only and is not intended for use in humans or animals.