

Rat Anti-mouse CD86/PE-Cy5 antibody

Catalog No: #28045



Package Size: #28045-1 50T #28045-2 100T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Anti-mouse CD86/PE-Cy5 antibody
Host Species	Rat
Clonality	Monoclonal
Applications	FC
Species Reactivity	Mouse
Specificity	Affinity purified by Protein G
Other Names	CD28LG2; FUN1; LAB72; Activation B7 2 antigen; Activation B7-2 antigen 3; CD86_HUMAN; CD86_MOUSE; Activation B7-2 antigen; Activation B72 antigen; B lymphocyte activation antigen B7 2; B lymphocyte activation antigen B72; B-lymphocyte activation antigen B7-2 2; B-lymphocyte activation antigen B7-2; B7 2; B7 2 antigen; B7; B7-2; B7.2; B70; B72 antigen; BU63; CD 86; CD28 antigen ligand 2 2; CD28 antigen ligand 2; Cd28l2; CD28LG2; CD86 antigen (CD28 antigen ligand 2, B7-2 antigen) 1, 2; CD86 antigen; CD86 molecule; CLS1; CTLA 4 counter receptor B7 2; CTLA 4 counter receptor B7.2; CTLA-4 counter-receptor B7.2 2, 3; CTLA4 counter receptor B72; Early T-cell costimulatory molecule 1; ETC-1; FUN 1; LAB72; Ly-58; Ly58; MB7; MB7-2; Membrane glycoprotein; MGC34413; T lymphocyte activation antigen CD86; T-lymphocyte activation antigen CD86; TS/A-2.
Accession No.	Swiss-Prot#: P42082 NCBI Gene ID:12524
Formulation	0.01M PBS, 0.5%BSA, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.

Application Details

Flow-Cyt=10ul/Test

Note: This product is for in vitro research use only and is not intended for use in humans or animals.