

PPM1M Polyclonal Antibody

Catalog No: #28580



Package Size: #28580-1 50ul #28580-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PPM1M Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse
Immunogen Description	Recombinant fusion protein of human PPM1M (NP_653242.3).
Conjugates	Unconjugated
Other Names	PPM1M; PP2C-eta; PP2CE; PP2Ceta; protein phosphatase 1M
Accession No.	Swiss-Prot#:Q96MI6NCBI Gene ID:132160
Calculated MW	30kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse testis, using PPM1M antibody.

Background

Published Papers

el at., Exploring the immunological role and prognostic potential of PPM1M in pan-cancerInMedicineOn2023 Mar 24byRongruo Zeng?1?2,?Lulu Wang et al..PMID:?36961170, , (2023)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.