

PKNOX1 Polyclonal Antibody

Catalog No: #28905



Package Size: #28905-1 50ul #28905-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PKNOX1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse
Immunogen Description	Recombinant fusion protein of human PKNOX1 (NP_004562.2).
Conjugates	Unconjugated
Other Names	PKNOX1; PREP1; pkonx1c; homeobox protein PKNOX1
Accession No.	Swiss-Prot#:P55347NCBI Gene ID:5316
Calculated MW	48kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of 293T cells, using PKNOX1 antibody.

Background

Published Papers

el at., The ring finger protein RNF6 induces leukemia cell proliferation as a direct target of pre-B-cell leukemia homeobox 1.In J Biol Chem on 2016 Apr 29 by Xin Xu, Kunkun Han,et al..PMID: 26971355, , (2016)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.