

TSPAN17 Polyclonal Antibody

Catalog No: #28958

Package Size: #28958-1 50ul #28958-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TSPAN17 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human TSPAN17 (NP_001006617.2).
Conjugates	Unconjugated
Other Names	TSPAN17; FBX23; FBXO23; TM4SF17; tetraspanin-17
Accession No.	Swiss-Prot#:Q96FV3NCBI Gene ID:26262
Calculated MW	35kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:200 - 1:2000

Images



Western blot analysis of extracts of mouse stomach, using TSPAN17 antibody.

Background

This gene encodes a member of the transmembrane 4 superfamily. It is characterized by four tetraspanin transmembrane segments. The function of this gene has not yet been determined.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.