

WDR82 Polyclonal Antibody

Catalog No: #29317



Package Size: #29317-1 50ul #29317-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	WDR82 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human WDR82 (NP_079498.2).
Conjugates	Unconjugated
Other Names	WDR82; MST107; MSTP107; PRO2730; PRO34047; SWD2; TMEM113; WDR82A; WD repeat domain 82
Accession No.	Swiss-Prot#:Q6UXN9NCBI Gene ID:80335
Calculated MW	35kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using WDR82 antibody.

Background

Published Papers

el at., Identification of Bone Morphometric Protein-Related Hub Genes and Construction of a Transcriptional Regulatory Network in Patients With Ossification of the Ligamentum FlavumInSpine (Phila Pa 1976)On2023 Sep 15byYifan Tuo?1?2,?Lihong Hu et al..PMID:?37384872, , (2023)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.