

MPHOSPH6 Polyclonal Antibody

Catalog No: #29416



Package Size: #29416-1 50ul #29416-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MPHOSPH6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human MPHOSPH6 (NP_005783.2).
Conjugates	Unconjugated
Other Names	MPHOSPH6; MPP; MPP-6; MPP6; M-phase phosphoprotein 6
Accession No.	Swiss-Prot#:Q99547NCBI Gene ID:10200
Calculated MW	22kDa,23kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MPHOSPH6 antibody.

Background

Published Papers

et al., 5-(3,4-Difluorophenyl)-3-(6-methylpyridin-3-yl)-1,2,4-oxadiazole (DDO-7263), a novel Nrf2 activator targeting brain tissue, protects against MPTP-Induced subacute Parkinson's disease in mice by inhibiting the NLRP3 inflammasome and protects PC12 cells against oxidative stress. In

Note: This product is for in vitro research use only and is not intended for use in humans or animals.