

Defb26 Polyclonal Antibody

Catalog No: #29483



Package Size: #29483-1 50ul #29483-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Defb26 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of rat Defb26 (NP_001032595.1).
Conjugates	Unconjugated
Other Names	DEFB126; C20orf8; DEFB-26; DEFB26; HBD26; bA530N10.1; hBD-26; beta-defensin 126
Accession No.	Swiss-Prot#:Q32ZG6NCBI Gene ID:641629
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IF 1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using Defb26 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.