

DHX57 Polyclonal Antibody

Catalog No: #29895



Package Size: #29895-1 50ul #29895-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DHX57 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Rat
Immunogen Description	Recombinant fusion protein of human DHX57 (NP_945314.1).
Conjugates	Unconjugated
Other Names	DHX57; DDX57; DExH-box helicase 57
Accession No.	Swiss-Prot#:Q6P158NCBI Gene ID:90957
Calculated MW	150-155kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DHX57 antibody.

Background

Published Papers

el at., High Levels of DEAH-Box Helicases Relate to Poor Prognosis and Reduction of DHX9 Improves Radiosensitivity of Hepatocellular Carcinoma. In Front Oncol on 2022 Jun 22 by Xi Chen, Letao Lin, et al..PMID: 35814441, , (2022)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.