

## GAPVD1 Polyclonal Antibody

Catalog No: #29981

Package Size: #29981-1 50ul #29981-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GAPVD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human GAPVD1 (NP_056450.2).
Conjugates	Unconjugated
Other Names	GAPVD1; GAPEX5; GAPex-5; RAP6; GTPase activating protein and VPS9 domains 1
Accession No.	Swiss-Prot#:Q14C86NCBI Gene ID:26130
Calculated MW	165kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:100

## Images



Western blot analysis of extracts of various cell lines, using GAPVD1 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.