

ABLIM2 Polyclonal Antibody

Catalog No: #30180

Package Size: #30180-1 50ul #30180-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABLIM2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse;Rat
Immunogen Description	Recombinant fusion protein of human ABLIM2 (NP_115808.3).
Conjugates	Unconjugated
Accession No.	Swiss-Prot#:Q6H8Q1NCBI Gene ID:84448
Calculated MW	68kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of Rat lung, using ABLIM2 antibody at 1:1000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.