

TRIM56 Polyclonal Antibody

Catalog No: #30308

Package Size: #30308-1 50ul #30308-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TRIM56 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human TRIM56.
Conjugates	Unconjugated
Accession No.	Swiss-Prot#:Q9BRZ2NCBI Gene ID:81844
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Background

Published Papers

Qiongzheng Zeng;Zixin Chen;Shan Li;Ziwei Huang;Zhe Ren;Cuifang Ye;Xiao Wang;Jun Zhou;Kaisheng Liu;Kai Zheng;Yifei Wang et al., HDAC6 deacetylates TRIM56 to negatively regulate cGAS-STING-mediated type I interferon responses., , (2025)

PMID:39747662

Note: This product is for in vitro research use only and is not intended for use in humans or animals.