

NEDD1 Polyclonal Antibody

Catalog No: #30637



Package Size: #30637-1 50ul #30637-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NEDD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human NEDD1 (NP_690869.1).
Conjugates	Unconjugated
Other Names	NEDD1; GCP-WD; TUBGCP7; protein NEDD1
Accession No.	Swiss-Prot#:Q8NHV4NCBI Gene ID:121441
Calculated MW	78kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using NEDD1 antibody at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using NEDD1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Published Papers

el at., NEDD1 overexpression increases cell proliferation, tumor immune escape, and drug resistance in LUAD. In J Cancer on 2024 Mar 11 by Ting Zhuo, Zuotao Wu,et al..PMID:38577589, , (2024)

[PMID:38577589](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.