

FUBP3 Polyclonal Antibody

Catalog No: #31385

Package Size: #31385-1 50ul #31385-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FUBP3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse
Immunogen Description	Recombinant fusion protein of human FUBP3 (NP_003925.1).
Conjugates	Unconjugated
Other Names	FUBP3;FBP3
Accession No.	Uniprot:Q96I24GeneID:8939
Calculated MW	62kDa
SDS-PAGE MW	62kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using FUBP3 antibody.

Background

May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.