

IL25 Polyclonal Antibody

Catalog No: #31486



Package Size: #31486-1 50ul #31486-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	IL25 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human IL25 (NP_073626.1).
Conjugates	Unconjugated
Other Names	IL25; IL17E; interleukin-25;IL-25
Accession No.	Swiss-Prot#:Q9H293NCBI Gene ID:64806
Calculated MW	24kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:2000
IHC□1:50 - 1:100
IF□1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using IL25 antibody.




Immunofluorescence analysis of human lung cancer cells using IL25 antibody.




Immunofluorescence analysis of rat bronchus cells using IL25 antibody.



Immunofluorescence analysis of mouse lung cells using IL25 antibody.



Immunohistochemistry of paraffin-embedded rat lung using IL25 antibody.



Immunohistochemistry of paraffin-embedded human prostate using IL25 antibody.

Background

The protein encoded by this gene is a cytokine that shares sequence similarity with interleukin 17. This cytokine can induce NF-kappaB activation, and stimulate the production of interleukin 8. Both this cytokine and interleukin 17B are ligands for the cytokine receptor IL17BR. Studies of a similar gene in mice suggest that this cytokine may be a pro-inflammatory cytokine favoring the Th2-type immune response. Alternative splicing results in multiple transcript variants.

Published Papers

el at., Interleukin-25 enhances humoral immune responses caused by the rabies virus. In Virulence on 2022 Dec by Yue Zhang, Mengwei Zhang,et al..PMID:35999776, , (2022)

[PMID:35999776](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.