

## MYBL1 Polyclonal Antibody

Catalog No: #31748

Package Size: #31748-1 50ul #31748-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MYBL1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human MYBL1 (NP_001138227.1).
Conjugates	Unconjugated
Other Names	MYBL1; A-MYB; AMYB; myb-related protein A
Accession No.	Swiss-Prot#:P10243NCBI Gene ID:4603
Calculated MW	110kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using MYBL1 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.