

## TTLL5 antibody

Catalog No: #31996

Package Size: #31996-1 50ul #31996-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TTLL5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Human;Mouse
Immunogen Description	Fusion protein of human TTLL5
Conjugates	Unconjugated
Target Name	TTLL5
Other Names	STAMP; CORD19; KIAA0998
Accession No.	NCBI Protein#:BC002766
Concentration	1.2mg/ml
Formulation	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Storage	Store at -20°C/1 year

## Application Details

IHC dilution:1:50-1:100

## Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using (TTLL5 Antibody) at dilution 1/70.

## Background

This gene encodes a member of the tubulin tyrosine ligase like protein family. This protein interacts with two glucocorticoid receptor coactivators, transcriptional intermediary factor 2 and steroid receptor coactivator 1. This protein may function as a coregulator of glucocorticoid receptor mediated gene induction and repression. This protein may also function as an alpha tubulin polyglutamylase.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.