

LAMP1 Antibody

Catalog No: #32731



Package Size: #32731-1 50ul #32731-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LAMP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of total LAMP1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human LAMP1 (NP_005552.3).
Conjugates	Unconjugated
Target Name	LAMP1
Other Names	LAMP1;CD107a;LAMPa;LGP120
Accession No.	Uniprot:P11279GeneID:3916
SDS-PAGE MW	90-120kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using LAMP1 antibody.



Immunofluorescence analysis of NIH-3T3 cells using LAMP1 Rabbit pAb.

Background

The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.

Published Papers

el at., Alamandine/MrgD axis prevents TGF- β 1-mediated fibroblast activation via regulation of aerobic glycolysis and mitophagyInJ Transl MedOn2023 Jan 13byWei Wang,Yue Zhang et al..PMID:36635651, , (2023)

[PMID:36635651](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.