

Anti-CD8a Antibody

Catalog No: #33199



Package Size: #33199-1 50ul #33199-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-CD8a Antibody
Host Species	Rat
Clonality	Monoclonal
Purification	The antibody was purified by affinity chromatography.
Applications	IF, FC
Species Reactivity	Mouse
Immunogen Description	Mouse thymus or spleen
Conjugates	Unconjugated
Other Names	MAL; T-cell surface glycoprotein CD8 alpha chain; T-lymphocyte differentiation antigen T8/Leu-2; CD8a
Accession No.	Swiss Prot:P01731GeneID:12525
Formulation	Rat IgG2a kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.