

MYB-A Antibody

Catalog No: #33667

Package Size: #33667-1 50ul #33667-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MYB-A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	The antibody detects endogenous levels of total MYB-A protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human MYB-A.
Conjugates	Unconjugated
Target Name	MYB-A
Other Names	AMYB; MYBA; MYBL1; Myb-related protein A;
Accession No.	Swiss-Prot: P10243NCBI Gene ID: 4603
SDS-PAGE MW	85kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from LOVO cells, using MYB-A antibody #33667.

Background

Strong transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Could have a role in the proliferation and/or differentiation of neurogenic, spermatogenic and B-lymphoid cells.

Nomura N., Nucleic Acids Res. 16:11075-11089(1988).

Golay J., Oncogene 9:2469-2479(1994).

Ma X.P., Cancer Res. 54:6512-6516(1994).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.