

NR4A3 Antibody

Catalog No: #34884



Package Size: #34884-1 50ul #34884-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NR4A3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous levels of total NR4A3 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human NR4A3.
Conjugates	Unconjugated
Target Name	NR4A3
Other Names	Nuclear receptor subfamily 4 group A member 3; Nuclear hormone receptor NOR-1; Neuron-derived orphan receptor 1; Mitogen-induced nuclear orphan receptor;
Accession No.	Swiss-Prot: Q92570NCBI Gene ID: 8013
SDS-PAGE MW	65kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from COLO cells, using NR4A3 antibody #34884.

Background

Binds to the B1A response-element.

Ohkura N., Biochim. Biophys. Acta 1308:205-214(1996).

Hedvat C.V., Mol. Endocrinol. 9:1692-1700(1995).

Clark J., Oncogene 12:229-235(1996).

Published Papers

el at., miR-665 expression predicts poor survival and promotes tumor metastasis by targeting NR4A3 In breast cancer. In Cell Death Dis. on 2019 Jun 17 by Zhao XG, Hu JY, et al..PMID:31209222, , (2019)

[PMID:31209222](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.