

## ECFP Mouse Polyclonal Antibody

Catalog No: #37991



Package Size: #37991-1 50ul #37991-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|               |                                                                                                                                                   |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name  | ECFP Mouse Polyclonal Antibody                                                                                                                    |
| Host Species  | Mouse                                                                                                                                             |
| Clonality     | Polyclonal                                                                                                                                        |
| Purification  | Affinity purification using immunogen.                                                                                                            |
| Applications  | WB                                                                                                                                                |
| Specificity   | The ECFP Antibody detects ECFP & CFP tag protein..                                                                                                |
| Target Name   | ECFP                                                                                                                                              |
| Concentration | 1.0mg/ml                                                                                                                                          |
| Formulation   | Mouse IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage       | Store at -20°C                                                                                                                                    |

## Application Details

Western blotting: 1:2000~1:5000

## Images



Western blot analysis of ECFP Recombinant protein, using #37991 diluted at 1) 1:5,000 2) 1:10,000

## Background

Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.