

β-actin Mouse Monoclonal Antibody

Catalog No: #38063



Package Size: #38063-1 50ul #38063-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	β-actin Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification using immunogen.
Applications	WB IHC IF
Species Reactivity	Zebrafish
Specificity	Antibody detects Zebrafish endogenous β-actin protein.
Conjugates	Unconjugated
Target Name	β-actin
Other Names	ACTB; BRWS1; PS1TP5BP1
Accession No.	Swiss-Prot#:P60709
SDS-PAGE MW	45kd
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:2000~1:5000

Images



Western blot analysis of Zebrafish skeletal muscle, using #38063 diluted at 1:3,000.

Background

β-actin is one of six different actin isoforms that have been identified. The actin molecules found in cells of various species and tissues tend to be very similar in their immunological and physical properties. Therefore, antibodies against β-actin are useful as loading controls for Western Blotting.

However it should be noted that levels of β-actin may not be stable in certain cells. For example, expression of β-actin in adipose tissue is very low

and therefore β -actin should not be used as loading control for these tissues.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.