

ECFP Rabbit Polyclonal Antibody

Catalog No: #38077

Package Size: #38077-1 50ul #38077-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ECFP Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB IP
Specificity	The ECFP Antibody detects ECFP & CFP tag protein..
Target Name	ECFP
Other Names	Cyan fluorescent protein; Enhanced Cyan fluorescent protein
Concentration	1.0mg/ml
Formulation	IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:2000~1:5000

Images



Western blot analysis of ECFP recombinant protein, using #38077 diluted at 1:5,000.

Background

Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.