

mOrange Rabbit Polyclonal Antibody

Catalog No: #38079

Package Size: #38079-1 50ul #38079-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	mOrange Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB IHC IF IP
Specificity	The mOrange Rabbit Polyclonal Antibody detects mOrange & mOrange Tag protein.
Target Name	mOrange
Other Names	mOrange
Concentration	1.0mg/ml
Formulation	IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:2000~1:5000

Images



Western blot analysis of mOrange recombinant protein, using #38079 diluted at 1:5,000.

Background

mOrange and mOrange2 are extremely bright orange fluorescent protein monomers which can be used as tags or reporters. Both mOrange fluorescent proteins are mutants derived from mRFP1, a monomeric mutant of DsRed.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.