

EBFP Mouse Polyclonal Antibody

Catalog No: #38081

Package Size: #38081-1 50ul #38081-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EBFP Mouse Polyclonal Antibody
Host Species	Mouse
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB IHC IF IP
Specificity	The EBFP Mouse Polyclonal Antibody detects EBFP & EBFP Tag protein.
Target Name	EBFP
Other Names	BFP; EBFP; Blue Fluorescent Protein; Enhanced Blue Fluorescent Protein
Accession No.	Swiss-Prot#:Q9ULX5
Concentration	1.0mg/ml
Formulation	Mouse IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:2000~1:5000

Images



Western blot analysis of EBFP recombinant protein, using #38081 diluted at 1:5,000.

Background

Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue fluorescence. EBFP can be used as protein “tags” to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.