

mStrawberry Mouse Monoclonal Antibody

Catalog No: #38099



Package Size: #38099-1 50ul #38099-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	mStrawberry Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification using immunogen.
Applications	WB
Specificity	The mStrawberry Mouse Polyclonal Antibody detects mStrawberry & mStrawberry Tag protein. It May had cross reactivity with mRFP protein.
Target Name	mStrawberry
Other Names	mStrawberry antibody
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:2000~1:5000

Images



Western blot analysis of mStrawberry recombinant protein, using #38099 diluted at 1:5,000.

Background

mStrawberry is a bright, monomeric red fluorescent protein which was developed in Dr. Roger Tsien's lab by directed mutagenesis of mRFP.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.