

DCX antibody

Catalog No: #38200

Package Size: #38200-1 50ul #38200-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DCX antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IF
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of total DCX protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human DCX.
Conjugates	Unconjugated
Target Name	DCX
Other Names	DCX;DBCN;DC;LISX;SCLH;XLIS;
Accession No.	Swiss-Prot#: O43602NCBI Gene ID: 1641
SDS-PAGE MW	45kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of mouse brain, using DCX antibody.

Background

Mutations in Doublecortin cause Lissencephaly (smooth brain), a neuronal migration disorder characterized by epilepsy and mental retardation (1).

Doublecortin is a microtubule associated protein that stabilizes and bundles microtubules. A conserved doublecortin domain mediates the interaction with microtubules, and interestingly most missense mutations cluster in this domain (2). Kinases JNK, CDK5 and PKA phosphorylate doublecortin. JNK phosphorylates Thr321, Thr331 and Ser334 while PKA phosphorylates Ser47 and CDK5 phosphorylates Ser297 (3-5). Phosphorylation of Ser297 lowers the affinity of doublecortin to microtubules. Furthermore, mutations of Ser297 result in migration defects (5).

Published Papers

el at., Naoxingqing tablet protects against cerebral ischemic/reperfusion injury by regulating ampka/NAMPT/SIRT1/PGC-1α pathwayInJ EthnopharmacolOn2024 Mar 25byXiao Sun?1,?Yunfeng Pan et al..PMID:?38159826, , (2023)

PMID:38159826

el at., Naoxingqing tablet protects against cerebral ischemic/reperfusion injury by regulating ampka/NAMPT/SIRT1/PGC-1α pathway. In J Ethnopharmacol on 2024 Mar 25 by Xiao Sun, Yunfeng Pan,et al..PMID:38159826, , (2024)

PMID:38159826

Note: This product is for in vitro research use only and is not intended for use in humans or animals.