

POMGNT2 antibody

Catalog No: #38956



Package Size: #38956-1 50ul #38956-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	POMGNT2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	IF
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total POMGNT2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human POMGNT2.
Target Name	POMGNT2
Other Names	AGO61; GTDC2; MDDGA8; C3orf39;
Accession No.	Swiss-Prot#: Q8NAT1NCBI Gene ID: 84892
SDS-PAGE MW	66kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:2000

Immunofluorescence: □1:50 - 1:100

Images



Immunofluorescence analysis of A549 cells using POMGNT2 antibody. Blue: DAPI for nuclear staining.

Background

This gene encodes a protein with glycosyltransferase activity although its function is not currently known.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.