

PRPF19 antibody

Catalog No: #39117

Package Size: #39117-1 50ul #39117-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PRPF19 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IHC
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of total PRPF19 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human PRPF19.
Conjugates	Unconjugated
Target Name	PRPF19
Other Names	PSO4; SNEV; PRP19; UBOX4; hPSO4; NMP200;
Accession No.	Swiss-Prot#: Q9UMS4NCBI Gene ID: 27339
SDS-PAGE MW	55kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:2000

Images



Western blot analysis of extracts of SW620 cells, using PRPF19 antibody at 1:1000 dilution.

Background

PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008)

Published Papers

el at., PRPF19 facilitates colorectal cancer liver metastasis through activation of the Src-YAP1 pathway via K63-linked ubiquitination of MYL9, , (2022)

[PMID:](#)

el at., PRPF19 facilitates colorectal cancer liver metastasis through activation of the Src-YAP1 pathway via K63-linked ubiquitination of MYL9InCell

Death DisOn2023 Apr 8byRui Zhou?#?1,?Jie Chen et al..PMID:?37031206, , (2023)

[PMID:37031206](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.