

Cyclin D1 Antibody

Catalog No: #39315

Package Size: #39315 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Cyclin D1 Antibody
Host Species	rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal Cyclin D1 antibody was raised against a 12 amino acid peptide from the C terminal residues of human Cyclin D1 protein.
Conjugates	Unconjugated
Target Name	Cyclin D1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded brain showing nuclear staining with anti-Cyclin D1 antibody dilution: 1:100.



Published Papers

Ye Liuqi;Lin Danlei;Zhang Wen;Chen Shiji;Zhen Yumiao;Akkouche Sara;Liang Xiaoxu;Chong Cheong-Meng;Zhong Hai-Jing; et al., AMBRA1 drives gastric cancer progression through regulation of tumor plasticity, , (2024)

[PMID:](#)

Zhipeng Zhang;Yanmei Wu;Wenqiang Liang;Zhifang Liao;Hongbo Liao;Xingxing Xing;Wenxin Yi;Zixuan Liu;Yicheng Li;Mengya Shi;Dongling Lin;Ting Gu;Biao Wu;Mingzhi Zou;Huilai Miao;Xin Wu et al., Eurycomalactone switched hepatocellular carcinoma cells into quiescence through 5'tRFAIa/DVL/ β -catenin pathway inhibition., , (2025)

[PMID:40128187](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.