

Osteocalcin Antibody

Catalog No: #39376



Package Size: #39376 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Osteocalcin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA IHC
Species Reactivity	Human;Mouse
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal osteocalcin antibody was raised against a peptide near the intermediate residues of human osteocalcin protein.
Conjugates	Unconjugated
Target Name	Osteocalcin
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA IHC(Paraffin section tested. Frozen section not been tested)

Published Papers

el at., Surface mineralized biphasic calcium phosphate ceramics loaded with urine-derived stem cells are effective in bone regeneration. In J Orthop Surg Res on 2019 Dec 9 by Xing , Li L, et al..PMID:31818319, , (2019)

[PMID:31818319](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.