

TNF α Antibody

Catalog No: #39389

Package Size: #39389 100ul



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TNF α Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB
Species Reactivity	Human
Specificity	Human
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein
Conjugates	Unconjugated
Target Name	TNF δO
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

WB

Images



Published Papers

Wang Ziheng;Lu Zhichao;Chen Yixun;Wang Chenxing;Gong Peipei;Jiang Rui;Liu Qianqian et al., Targeting the AKT-P53/CREB pathway with epicatechin for improved prognosis of traumatic brain injury, , (2023)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.