

XCR1 Antibody

Catalog No: #39712

Package Size: #39712 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	XCR1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Immunogen Description	raised against a peptide from C terminal residues of human XCR1 protein.
Target Name	XCR1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA WB IHC

Images



Published Papers

el at., Transcriptomic analyses of chemokines reveal that down-regulation of XCR1 is associated with advanced hepatocellular carcinoma .On Biochem Biophys Res Commun. 2018 Feb 19 by Yanru W, Zhenyu B et al..PMID: 29408492, , (2018)

[PMID:29408492](https://pubmed.ncbi.nlm.nih.gov/29408492/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.