

VNN2 Antibody

Catalog No: #39964

Package Size: #39964 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VNN2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal VNN2(1) antibody was raised against a recombinant human VNN2 protein 254-435aa (BC126145).
Conjugates	Unconjugated
Target Name	VNN2
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% sodium azide.
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA: 1:20000-1:80000 WB: 1:200-1:1000 IHC: 1:100-1:500

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded breast carcinoma showing membrane staining with anti-VNN2 antibody at a dilution of 1/100



Published Papers

et al., Upregulated Vanins and their potential contribution to periodontitis. In BMC Oral Health on 2022 Dec 17 by Weijun Yu, Shucheng Hu, et al..PMID: 36527111, , (2022)

[PMID:36527111](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.