

MOCS3 Antibody

Catalog No: #40000

Package Size: #40000 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MOCS3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal MOCS3 (1) antibody was raised against a recombinant human MOCS3 protein 15-175aa (BC015939).
Conjugates	Unconjugated
Target Name	MOCS3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% sodium azide.
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

WB 1:500-2000

Images



Published Papers

Mi Yuan, Mao Junjie, Li Xuzhe, Wang Zilin, Zuo Xin, Chang Zihan, Wang Ruibo, Wang Lantao, Wang Lei et al., AURKA modulates tight junction functionality to influence the proliferation and metastasis of lung adenocarcinoma, *Molecular and cellular biochemistry*, (2025)

[PMID:41201786](https://pubmed.ncbi.nlm.nih.gov/41201786/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.