

GFPT1 Rabbit mAb

Catalog No: #40024

Package Size: #40024 100ul



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GFPT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Type	peptide
Immunogen Description	A synthesized peptide derived from human GFPT1
Conjugates	Unconjugated
Target Name	GFPT1
Concentration	1.0mg/ml
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Application Details

WB 1:1000 - 1:5000;IHC 1:100 - 1:200

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.